

Dantec™ Keypoint® Focus

Compact EMG/NCS/EP Workstation

The Dantec Keypoint Focus EMG system is based on nearly 50 years of experience in EMG system design. Quality, accuracy, and dependability – what you have come to expect from the Dantec brand.



Dantec Keypoint Focus is the next generation EMG system – delivering the quality and results you expect, at an economical price.

Its new compact design and intuitive, powerful functionality make the Dantec Keypoint Focus EMG system ideal for space conscious professionals in hospitals and private offices.

- All-in-One EMG/EP system with dedicated controls, integrated stimulators and high quality audio
- 3-, 4-, 6- or 8-channel amplifier options
- Advanced features for quantitative EMG, SF-EMG and tremor analysis
- On-line comparison to reference values

Handle difficult cases more quickly and efficiently

- Protocol-based testing
- Broad range of powerful analysis tools
- Automatic comparison with normative data
- The only system available with a complete MUP normative database
- Fast, automatic report generation

Design an EMG system that's perfect for you

- Select the configuration that works for you - portable notebook systems to fully featured workstations
- Purchase only the applications you need
- Completely control all measurements, data tables/calculations and displays
- Create customized reports


IT-friendly features

- Microsoft SQL database
- Automated study data back-up/archiving functionality
- Complete Active Directory Services
- GDT/HL7 connectivity

Versatile EMG/NCS/EP Software




Multimodality NCS with reference values




Advanced H Reflex




Advanced Decrement test



Multi-MUP with reference values



Peak-triggering SF-EMG



EP with reference values

Natus Neuro Support

Excellence in customer service and support is a core value of Natus Neuro.

Essential Support Elements:

- Accessible and effective Technical Support
- Definitive technical documentation and knowledgeable installation teams
- Excellent replacement unit or spare part availability
- Extended warranty and service coverage programs
- Comprehensive, flexible customer training courses

Supplies

Natus Neuro offers a full range of neurodiagnostic supplies.

- Dedicated and knowledgeable customer service
- Streamlined order processing to save you time



Natus Europe GmbH

Robert-Koch-Str. 1
82152 Planegg
Germany
Tel: +49 (0)89 / 83942-0
Fax: +49 (0)89 / 83942-777

neuro.natus.com

Natus Neuro

3150 Pleasant View Road
Middleton, WI 53562 USA

To learn more about Natus Neuro products contact your local distributor or sales representative.

Call: **+1-608-829-8500**