



Dantec Clavis

Handheld device for EMG/STIM-guided injections

Simple sophistication

Clavis is designed to be simple and intuitive to provide accurate guidance in minimal time

A fast & easy pathway to injection of medicinal products

Clavis is designed to provide the functionality of EMG recording and stimulation with mobility and simplicity. Clavis combines the quality of larger systems with the ease and flexibility of a handheld device.

- EMG and STIM functions in one compact device
- Easy to use and operate
- Battery powered
- Perfect match with Bo-ject DHN needle electrodes

Disposable Bo-ject® DHN disposable hypodermic needle electrodes

| Part Number | Needle Length | Needle Diameter | Lead Color | Quantity |
|-------------|---------------|-----------------|-------------|------------------|
| 901350422 | 1" (25mm) | 0.30mm (30G) | Tan | 1/pouch - 10/box |
| 901350432 | 1" (25mm) | 0.41mm (27G) | Pink | 1/pouch - 10/box |
| 901350442 | 1.5" (37mm) | 0.41mm (27G) | Light Blue | 1/pouch - 10/box |
| 901350472 | 1.5" (37mm) | 0.46mm (26G) | Brown | 1/pouch - 10/box |
| 901350452 | 2" (50mm) | 0.51mm (25G) | Grey | 1/pouch - 10/box |
| 901350462 | 3" (75mm) | 0.71mm (22G) | Light Green | 1/pouch - 10/box |

Disposable pre-gelled surface electrodes


| Part # | Recording Area | Lead Length | Lead Connector | Quantity |
|--------|----------------|---------------|-------------------------|----------|
| 443010 | 15mm x 20mm | 6.5" (16.5cm) | 0.7mm male touchproof | 12/pkg |
| 454010 | 20mm x 25mm | 6.5" (16.5cm) | 0.7mm male touchproof | 12/pkg |
| 443020 | 15mm x 20mm | 47.2" (1.2m) | 1.5mm female touchproof | 12/pkg |
| 454030 | 20mm x 25mm | 47.2" (1.2m) | 1.5mm female touchproof | 12/pkg |
| 454040 | 20mm x 25mm | 79" (2.0m) | 1.5mm female touchproof | 12/pkg |

Disposable ground/tab electrode

| Part # | Type | Quantity |
|------------|--------|----------------|
| 019-435300 | Ground | 10 pouches/box |

Alligator clip cable for tab electrode

| Part # | Type | Length | Electrode End | Instrument End | Quantity |
|-----------|------------|------------|---------------------------|-------------------------------|----------|
| 9013P0602 | Unshielded | 31" (80cm) | Clip Connector – Green | 1.5mm touchproof connector | 1/pkg |
| 9013P0642 | Unshielded | 31" (80cm) | Clip Connector – Red | 1.5mm touchproof connector | 1/pkg |



Benefits of a large system in a handheld device

Clavis provides muscle and nerve localization wherever needed to accurately guide and monitor injection of medicinal products. It provides EMG recording with audio and stimulation ranging from 0-15mA with speed and reliability while efficiently operating on battery power.

Guidance for injection of medicinal products where & when you need it

Clavis is quick and easy to operate. Provides guidance for accurate injection of medicinal products at the exact point of treatment.

Clavis is small and lightweight for easy, portable use in any location.

Instantly switch between stimulation & recording

Switch quickly and easily between stimulation and recording functions.

When seconds really count, Clavis' fast and simple operation saves time.



Clavis & Bo-ject needles – made for each other

Bo-ject needle electrodes are:

- Consistently sharp, stainless steel needle with a low-friction coating for patient comfort
- Compatible with a wide range of reference electrodes and syringes with Luer Lock or slip tip connectors
- Easy to attach to Clavis and ready to use quickly
- Ergonomically designed and lightweight for greater control



To ensure consistently reliable results, Clavis is designed to work with Bo-ject disposable hypodermic needle electrodes.

These needle electrodes have been specifically developed for injection of medicinal products such as botulinum toxin.



Service

Excellence in customer and technical service is a core value of Natus Neurology.

Essential Support Elements:

- Accessible and effective Technical Support
- Definitive technical documentation and web-training
- Excellent replacement/exchange availability

Supplies

Natus Neurology offers a full range of neurodiagnostic supplies.

- Dedicated and knowledgeable customer service
- Streamlined order processing to save you time
- Convenient online portal (US Customers only)
[Natus Neuro Store – www.natusneurostore.com](http://www.natusneurostore.com)

To learn more about Natus Neurology Service Programs or our full line of Supplies, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**

natus
neurology

Natus Neurology Incorporated

3150 Pleasant View Road
Middleton, WI 53562 USA
Tel: 1-800-356-0007
1-608-829-8500
Fax: 1-608-829-8709
www.natus.com

©2017 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.