



## Genuine Grass<sup>®</sup> Reusable Cup Electrodes

### Finest Materials, Reliable Design, Precision Craftsmanship

- Every Genuine Grass Electrode is tested for mechanical and electrical properties to assure long life... backed by a ONE-YEAR\* warranty
- Unsurpassed recording clarity and dependability
- Choice of deep, 10mm diameter, cast cups in gold or silver or 6mm diameter gold cast cups
- Choice of three unique, versatile 1.5mm touchproof connectors
- Choice of lightweight or tangle-free lead wires; both include Kevlar<sup>®</sup> strands for added strength
- Compatible with steam and EtO sterilization

Heavy 24k gold plate over  
pure silver cast cup



Pure, solid silver cast cup,  
burnished for a smooth finish



Heavy 24k gold plate  
over pure silver



Hole drilled in stem permits full lead wire insertions for optimum electrical and mechanical bond. Double strain relieved for durability and to prevent fluid invasion to maintain purity.

# Genuine Grass Reusable Cup Electrodes

Item Number	Cup Type	Wire Length & Type	Connector Type	Quantity
F-E5GH-12	10mm Gold w/2mm hole	12" Lightweight Teflon	Molded 1.5mm Touchproof	10/box
F-E5GH-24	10mm Gold w/2mm hole	24" Lightweight Teflon	Molded 1.5mm Touchproof	10/box
F-E5GH-30	10mm Gold w/2mm hole	30" Lightweight Teflon	Molded 1.5mm Touchproof	10/box
F-E5GH-48	10mm Gold w/2mm hole	48" Lightweight Teflon	Molded 1.5mm Touchproof	10/box
F-E5GH-60	10mm Gold w/2mm hole	60" Lightweight Teflon	Molded 1.5mm Touchproof	10/box
F-E5GH-72	10mm Gold w/2mm hole	72" Lightweight Teflon	Molded 1.5mm Touchproof	10/box
F-E5GH-96	10mm Gold w/2mm hole	96" Lightweight Teflon	Molded 1.5mm Touchproof	10/box
F-E5GH-120	10mm Gold w/2mm hole	120" Lightweight Teflon	Molded 1.5mm Touchproof	10/box
F-E6GH-48	6mm Gold w/1.5mm hole	48" Lightweight Teflon	Molded 1.5mm Touchproof	10/box
F-E5SH-48	10mm Silver w/2mm hole	48" Lightweight Teflon	Molded 1.5mm Touchproof	10/box
F-E5SH-72	10mm Silver w/2mm hole	72" Lightweight Teflon	Molded 1.5mm Touchproof	10/box
FH-E5GH-48	10mm Gold w/2mm hole	48" Lightweight Teflon	Hand-Crafted 1.5mm Touchproof	10/box
FH-E5GH-60	10mm Gold w/2mm hole	60" Lightweight Teflon	Hand-Crafted 1.5mm Touchproof	10/box
FH-E5GH-72	10mm Gold w/2mm hole	72" Lightweight Teflon	Hand-Crafted 1.5mm Touchproof	10/box
FH-E5GH-96	10mm Gold w/2mm hole	96" Lightweight Teflon	Hand-Crafted 1.5mm Touchproof	10/box
FH-E5SH-48	10mm Silver w/2mm hole	48" Lightweight Teflon	Hand-Crafted 1.5mm Touchproof	10/box
FD-E5GH-48	10mm Gold w/2mm hole	48" Lightweight Teflon	Extra-Width 1.5mm Touchproof	10/box
FD-E5GH-60	10mm Gold w/2mm hole	60" Lightweight Teflon	Extra-Width 1.5mm Touchproof	10/box
FD-E5GH-72	10mm Gold w/2mm hole	72" Lightweight Teflon	Extra-Width 1.5mm Touchproof	10/box
FD-E5GH-96	10mm Gold w/2mm hole	96" Lightweight Teflon	Extra-Width 1.5mm Touchproof	10/box
FS-E5GH-48	10mm Gold w/2mm hole	48" No-Tangle Silicone	Molded 1.5mm Touchproof	10/box
FS-E5GH-60	10mm Gold w/2mm hole	60" No-Tangle Silicone	Molded 1.5mm Touchproof	10/box
FS-E5GH-72	10mm Gold w/2mm hole	72" No-Tangle Silicone	Molded 1.5mm Touchproof	10/box
FS-E5GH-96	10mm Gold w/2mm hole	96" No-Tangle Silicone	Molded 1.5mm Touchproof	10/box
FS-E5GH-120	10mm Gold w/2mm hole	120" No-Tangle Silicone	Molded 1.5mm Touchproof	10/box
FS-E5SH-48	10mm Silver w/2mm hole	48" No-Tangle Silicone	Molded 1.5mm Touchproof	10/box
FS-E5SH-60	10mm Silver w/2mm hole	60" No-Tangle Silicone	Molded 1.5mm Touchproof	10/box

**Ten Lead Wire Colors:** Brown, Red, Orange, Yellow, Green, Blue, Violet, Gray, White, Tan

## Materials

### Electrode

Gold Cup ..... Heavy 24k gold plate over pure silver, burnished for smoothness  
Silver Cup ..... Pure silver, burnished for smoothness

### Lead Wire

#### Conductor

Compliant, non-kink 7-strands of silver plated copper alloy with Kevlar core for strength

#### Insulation

Highly flexible Teflon (F, FH, FD Series) or extremely supple tangle-free Silicone (FS Series)  
Polyolefin heat-shrink tubing to prevent fluid invasion

### Touchproof connector, 1.5mm (DIN 42 802)

Connector..... Copper alloy body, stainless steel sleeve

Insulation..... Pelprene thermoplastic resin (F, FS Series); Polyolefin (FH, FD Series)

Latex Free / RoHS compliant



FH Series hand-crafted touchproof connector



FD Series extra-width touchproof connector



F and FS Series molded plastic touchproof connector

## Supplies

Natus Neurology offers a full range of neurodiagnostic supplies promoting patient comfort. Our dedicated customer service team provides a streamlined order and shipping process to save you time and money.

To learn more about Natus Neurology full line of Supplies, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**

**natus**  
neurology

#### Natus Neurology Incorporated

3150 Pleasant View Road  
Middleton, WI 53562 USA

Tel: 1-800-356-0007

1-608-829-8500

Fax: 1-608-829-8709

[www.natus.com](http://www.natus.com)

#### Natus Manufacturing Limited

IDA Business Park, Gort  
Co. Galway, Ireland.

© 2014 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.