

**New  
25-packs**

# Grass<sup>®</sup> Disposable Deep Cup Electrodes

(Silver/Silver Chloride, 10mm)

## Trusted, Reliable, Convenient

**Grass** Disposable Deep Cup Electrodes continue a tradition of quality and innovation from the most trusted brand in EEG electrodes. Increased cup volume holds more paste for extended EEG recording.

- 50% more surface contact area and 20% greater cup depth than standard models
- Low impedance with both routine EEG and long term monitoring
- No cleaning and maintenance; use once and dispose
- Touchproof connector fit is optimized for easy headbox insertion and removal
- Packs of 10: ten unique lead wire colors  
●●●●●●●●●●
- Packs of 25: twelve unique lead wire colors  
(two each plus one extra white) ●●●●●●●●●●●●

To learn more about Natus Neuro full line of Supplies, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**

Grass Disposable Deep Cup Electrodes			
Part Number	Cup Type	Lead Length	Quantity
023146	10mm Ag/AgCl	39" (1.0m)	10/pkg
023147	10mm Ag/AgCl	59" (1.5m)	10/pkg
023148	10mm Ag/AgCl	79" (2.0m)	10/pkg
023149	10mm Ag/AgCl	98" (2.5m)	10/pkg
027924	10mm Ag/AgCl	39" (1.0m)	25/pkg
027925	10mm Ag/AgCl	59" (1.5m)	25/pkg
027926	10mm Ag/AgCl	79" (2.0m)	25/pkg
027927	10mm Ag/AgCl	98" (2.5m)	25/pkg

## Materials

### Electrode

Silver/Silver Chloride coated ABS plastic

### Cup Dimensions

Cup diameter ..... 10mm  
Cup height ..... 3.0mm  
Interior cup diameter ..... 6.7mm

### Lead Wire

Core ..... Tin plated Copper  
Jacket ..... PVC

### Touchproof connector, 1.5mm female (DIN 42 802)

Connector ..... Tin plated Brass  
Insulation ..... PVC

**RoHS compliant**

New  
25-packs

# Grass<sup>®</sup> Disposable Standard Cup Electrodes

(Silver/Silver Chloride, 10mm)

## Trusted, Reliable, Convenient

**Grass** Disposable Standard Cup Electrodes continue a tradition of quality and innovation from the most trusted brand in EEG electrodes. The standard cup electrodes are designed with a wide, flat rim for gentle contact with sensitive skin.

- Low impedance with both routine EEG and long term monitoring
- No cleaning and maintenance; use once and dispose
- Touchproof connector fit is optimized for easy headbox insertion and removal
- Packs of 10: ten unique lead wire colors  
●●●●●●●●●●
- Packs of 25: twelve unique lead wire colors  
(two each plus one extra white) ●●●●●●●●●●●●

**natus**

Natus Medical Incorporated  
3150 Pleasant View Road  
Middleton, WI 53562 USA  
[natus.com](http://natus.com)

Grass Disposable Standard Cup Electrodes			
Part Number	Cup Type	Lead Length	Quantity
019-477200	10mm Ag/AgCl	39" (1.0m)	10/pkg
019-477300	10mm Ag/AgCl	59" (1.5m)	10/pkg
019-477400	10mm Ag/AgCl	79" (2.0m)	10/pkg
019-479200	10mm Ag/AgCl	98" (2.5m)	10/pkg
027920	10mm Ag/AgCl	39" (1.0m)	25/pkg
027921	10mm Ag/AgCl	59" (1.5m)	25/pkg
027922	10mm Ag/AgCl	79" (2.0m)	25/pkg
027923	10mm Ag/AgCl	98" (2.5m)	25/pkg

## Materials

### Electrode

Silver/Silver Chloride coated ABS plastic

### Cup Dimensions

Cup diameter ..... 10mm  
Cup height ..... 2.5mm  
Interior cup diameter ..... 4.6mm

### Lead Wire

Core ..... Tin plated Copper  
Jacket ..... PVC

### Touchproof connector, 1.5mm female (DIN 42 802)

Connector ..... Tin plated Brass  
Insulation ..... PVC

### RoHS compliant

To learn more about Natus Neuro full line of Supplies, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**