



Grass Reusable “Tangle Free” Cup Electrodes

Durability, Ease of Use and Reliable Signals Time After Time

- Choice of gold, silver or silver/silver chloride 10mm cups
- Pediatric 6mm cup silver/silver chloride available
- No-lift contour adjustable cup neck for low profile (gold & silver cups only)
- Molded plastic 1.5mm female touchproof connector
- Low insertion force connector for easy headbox insertion and removal
- Six lead wire colors (two each): red, blue, green, white, brown, yellow

Part Number	Cup Type	Lead Length	Quantity
019-477500	10mm Gold	39" (1.0m)	12/pkg
019-477600	10mm Gold	59" (1.5m)	12/pkg
019-478600	10mm Gold	79" (2.0m)	12/pkg
019-477700	10mm Gold	98" (2.5m)	12/pkg
019-477800	10mm Silver	39" (1.0m)	12/pkg
019-477900	10mm Silver	59" (1.5m)	12/pkg
019-478700	10mm Silver	79" (2.0m)	12/pkg
019-478800	10mm Silver	98" (2.5m)	12/pkg
019-478400	10mm Ag/AgCl	39" (1.0m)	12/pkg
019-478500	10mm Ag/AgCl	59" (1.5m)	12/pkg
019-478900	10mm Ag/AgCl	79" (2.0m)	12/pkg
019-479000	10mm Ag/AgCl	98" (2.5m)	12/pkg
019-479100	6mm Ag/AgCl Pediatric	59" (1.5m)	12/pkg

Materials

Electrode

Gold, Silver or Silver/Silver Chloride

Lead Wire

Core Kevlar and Silver plated Copper
Jacket TPU

Touchproof connector, 1.5mm female (DIN 42 802)

Connector Gold plated Brass
Insulation PVC

Latex Free / RoHS compliant

Grass Disposable Cup Electrodes (Silver/Silver Chloride, 10mm)

Safety, Efficiency and Reliable Signals

- Single use eliminates the risk of cross-contamination
- No cleaning and maintenance; use once and dispose
- Flexible lead wires terminate in 1.5mm female touchproof connectors
- Five lead wire colors (two each): green, blue, red, yellow, white

Part Number	Cup Type	Lead Length	Quantity
019-477200	10mm Ag/AgCl	39" (1.0m)	10/pkg
019-477300	10mm Ag/AgCl	59" (1.5m)	10/pkg
019-477400	10mm Ag/AgCl	79" (2.0m)	10/pkg
019-479200	10mm Ag/AgCl	98" (2.5m)	10/pkg

Materials

Electrode

Silver/Silver Chloride coated ABS plastic

Lead Wire

Core Tin plated Copper


Jacket PVC

Touchproof connector, 1.5mm female (DIN 42 802)

Connector Tin plated Brass

Insulation PVC

Latex Free / RoHS compliant



Supplies

Natus Neurology offers a full range of neurodiagnostic accessories and supplies promoting patient comfort. Our dedicated customer service team provides a streamlined order and shipping process to save you time and money.

To learn more about Natus Neurology full line of Supplies, contact your local distributor or sales representative.

US Customers Call:
1-800-356-0007

International Customers Call:
+1-608-829-8500



Natus Neurology Incorporated

3150 Pleasant View Road
Middleton, WI 53562 USA

Tel: 1-800-356-0007

1-608-829-8500

Fax: 1-608-829-8709

www.natus.com

© 2015 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.