



## Natus Ultra Subdermal Needle Electrodes

Quality... Comfort... Value



- Market pricing
- Excellent quality; ease of insertion; reliable signal
- Trusted source for neurodiagnostic instruments and IOM supplies

- Precision cut, gently tapered, ultra-sharp lancet tip for effortless skin penetration assuring optimal patient comfort
- High-tensile-strength, medical grade stainless steel needle for optimum strength and flexibility
- Low profile needle-to-wire joint reduces pressure points on the patient's skin
- Multi-colored lead wires with matching connectors for easy placement identification
- Low noise and impedance level for clear, reliable recordings
- 1.5mm insulated female connector, touchproof (DIN 42 802)
- Individually packaged in sterile, pre-printed pouches for safety and easy identification



## Natus Ultra Subdermal Needle Electrodes

Item Number	Type	Needle Length/ Diameter	Lead Wire Length	Quantity
019-476600	Stainless Steel, 6-color	13mm, 27g(0.38mm)	1m/39"	1/pouch, 24/box
019-476800	Stainless Steel, 6-color	13mm, 27g(0.38mm)	1.5m/59"	1/pouch, 24/box
019-477000	Stainless Steel, 6-color	13mm, 27g(0.38mm)	2.5m/96"	1/pouch, 24/box
019-476700	Stainless Steel, 12-color	13mm, 27g(0.38mm)	1m/39"	1/pouch, 24/box
019-476900	Stainless Steel, 12-color	13mm, 27g(0.38mm)	1.5m/59"	1/pouch, 24/box
019-477100	Stainless Steel, 12-color	13mm, 27g(0.38mm)	2.5m/96"	1/pouch, 24/box

4 each of 6 Lead Wire Colors: Red, Green, White, Black, Blue, Yellow



2 each of 12 Lead Wire Colors: Red, Green, White, Black, Blue, Yellow, Gray, Orange, Brown, Purple, Light Blue, Pink



### Materials:

Needle ..... Medical grade Stainless Steel  
 Wire..... Copper  
 Wire Insulation.....PVC  
 Connector Type ..... Touchproof, 1.5mm insulated female (DIN 42 802)  
 Sterilization Type .....EtO  
 Shelf Life ..... 3 years for unopened pouches



Natus Ultra Subdermal touchproof connector

## Supplies

Natus Neurology offers a full range of neurodiagnostic accessories and supplies promoting patient comfort. Our dedicated customer service team provides a streamlined order and shipping process to save you time and money.

To learn more about the full line of Natus Neurology Supplies contact your local distributor or sales representative:

US/Canada Customers Call:  
**1-800-356-0007**

International Customers:  
**Contact your local distributor**  
View our full line of supplies at:  
**natus.com/supplies**



### Natus Neurology Incorporated

3150 Pleasant View Road  
 Middleton, WI 53562 USA  
 Tel: 1-800-356-0007  
 1-608-829-8500  
 Fax: 1-608-829-8709  
[www.natus.com](http://www.natus.com)

© 2015 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.