







Natus® XactTrace® Respiratory Effort Belts

RIP* effort belt solutions sized to fit a wide range of patients

The Natus XactTrace family of respiratory effort belts provide a wide range of sizing options for in-home and in-lab sleep studies. XactTrace belts are reliable and have improved cabled strength relief for increased durability.

- 1 **Reusable Belts** provide quality data, are economical, and fit most patients
- 2 **Single Use Cut-To-Fit Belts** allow customized sizing to accommodate hard-to-fit patients
- 3 **Single Use Pre-sized Belts** provide ease of use, fast application and quality data

1 Reusable Belts			
	XactTrace Reusable Universal RIP Belts	Compatible with Xltek, Grass, Nicolet, Sandman SD32 and most other PSG systems that accept 1.5mm female touchproof connectors	<ul style="list-style-type: none"> • Integrated battery power • Adjustable • Fits most patients
	XactTrace Reusable Embla RIP Belts	Compatible with Embla N7000, Embla S4000/S4500, Embla titanium and all Embletta systems	<ul style="list-style-type: none"> • Powered by Embla PSG System • Small, Medium and Large sizes
2 Single Use Cut-To-Fit Belts			
	XactTrace Single Use Cut-To-Fit Universal RIP Belts	Compatible with Xltek, Grass, Nicolet, Sandman SD32 and most other PSG systems that accept 1.5mm female touchproof connectors	<ul style="list-style-type: none"> • Powered by an external RIP processor • Customized sizing for each patient; ideal for pediatric and bariatric use • Use only Cut-To-Fit Locks/Snap Sensors
	XactTrace Single Use Cut-To-Fit Embla RIP Belts	Compatible with Embla N7000, Embla S4000/S4500, Embla titanium, Embletta MPR, Embletta Gold**, Embletta 128 and Embletta X100 systems	<ul style="list-style-type: none"> • Powered by Embla PSG System • Customized sizing for each patient; ideal for pediatric and bariatric use • Use only Cut-To-Fit Locks/Snap Sensors • Optional Cable Extensions
3 Single Use Pre-Sized Belts			
	XactTrace Single Use Pre-Sized Universal RIP Belts	Compatible with most PSG systems that accept 1.5mm female touchproof connectors	<ul style="list-style-type: none"> • Powered by an external RIP processor • Pediatric, Small, Medium and Large sizes • Use only Pre-Sized Locks/Snap Sensors
	XactTrace Single Use Pre-Sized Embla RIP Belts	Compatible with Embla N7000, Embla S4000/S4500, Embla titanium, Embletta MPR (all models) and Embletta Gold*** systems	<ul style="list-style-type: none"> • Powered by Embla PSG System • Pediatric, Small, Medium and Large sizes • Use only Pre-Sized Locks/Snap Sensors

*Respiratory Inductance Plethysmography **Embletta Gold requires locks 1421021 and 1421023 ***Embletta Gold requires 2 of part 013382

Ordering Information

Additional XactTrace product information is available in the **Neurodiagnostic Supplies catalog** and on the Natus website www.natus.com/supplies

1 XactTrace Reusable UNIVERSAL RIP Belt Kits

Part Number	Type/Description
1421091	Universal Kit - Complete
1421080	Universal Abdomen Kit (belt, sensor, cable)
1421081	Universal Thorax Kit (belt, sensor, cable)

XactTrace Reusable RIP Belts for EMBLA PSG Systems

Part Number	Type/Description
1421032	Small ABD, 22" (57cm) to 41" (104cm), 27" (70cm) cable
1421033	Small TRX, 22" (57cm) to 41" (104cm), 15.7" (40cm) cable
1421034	Medium ABD, 33" (83cm) to 56" (142cm), 27.5" (70cm) cable
1421035	Medium TRX, 33" (83cm) to 56" (142cm), 15.7" (40cm) cable
1421036	Large ABD, 47" (120cm) to 76" (194cm), 27.5" (70cm) cable
1421037	Large TRX, 47" (120cm) to 76" (194cm), 15.7" (40cm) cable
1421042	Small ABD, 22" (57cm) to 41" (104cm), 78" (200cm) cable
1421041	Small TRX, 22" (57cm) to 41" (104cm), 78" (200cm) cable
1421039	Medium ABD, 33" (83cm) to 56" (142cm), 78" (200cm) cable
1421038	Medium TRX, 33" (83cm) to 56" (142cm), 78" (200cm) cable
1421044	Large ABD, 47" (120cm) to 76" (194cm), 78" (200cm) cable
1421043	Large TRX, 47" (120cm) to 76" (194cm), 78" (200cm) cable

2 XactTrace Single Use Cut-To-Fit UNIVERSAL RIP Belt Kits

Part Number	Type/Description
1421052	Universal Single Use Cut-to-Fit Belt Kit - Complete*
1421050	Universal Single Use Cut-to-Fit Abdomen Kit*
1421051	Universal Single Use Cut-to-Fit Thorax Kit*

XactTrace Single Use Cut-To-Fit RIP Belts Components for EMBLA PSG Systems

Part Number	Type/Description
1421020	Lock for Cut-To-Fit Single Use Belt, ABD, 27.5" (70cm) cable
1421022	Lock for Cut-To-Fit Single Use Belt, Thorax, 15.7" (40cm) cable
1421021	Lock for Cut-To-Fit Single Use Belt, ABD, 27.5" (70cm) cable**
1421023	Lock for Cut-To-Fit Single Use Belt, Thorax, 15.7" (40cm) cable**
1421000	Cut-To-Fit Single Use belt material, 65.5' (20m)
1451070	Cable Extension, Abdomen, 78" (200cm), Optional
1451072	Cable Extension, Thorax, 78" (200cm), Optional

*Kits include locks, RIP processors and belt material

**For Embletta Gold only

3 XactTrace Single Use Pre-Sized UNIVERSAL RIP Belts, Snap Sensors & RIP Processors

Part Number	Type/Description
011516	Pre-Sized Single Use Belt, Pediatric
011517	Pre-Sized Single Use Belt, Small
011518	Pre-Sized Single Use Belt, Medium
011520	Pre-Sized Single Use Belt, Large
011521	Snap Sensor for Pre-Sized Single Use Belts, Abdomen, 26" (66cm) cable
011522	Snap Sensor for Pre-Sized Single Use Belts, Thorax, 16" (41cm) cable
1421048	RIP Processor Sensor, Abdomen, 78" (200cm) cable
1421047	RIP Processor Sensor, Thorax, 78" (200cm) cable

XactTrace Single Use Pre-Sized RIP Belts & Snap Sensors for EMBLA PSG Systems

Part Number	Type/Description
011516	Pre-Sized Single Use Belt, Pediatric
011517	Pre-Sized Single Use Belt, Small
011518	Pre-Sized Single Use Belt, Medium
011520	Pre-Sized Single Use Belt, Large
011521	Snap Sensor for Pre-Sized Single Use Belts, Abdomen, 26" (66cm) cable
011522	Snap Sensor for Pre-Sized Single Use Belts, Thorax, 16" (41cm) cable
013224	Snap Sensor for Pre-Sized Single Use Belts, Abdomen, 78" (200cm) cable
013225	Snap Sensor for Pre-Sized Single Use Belts, Thorax, 78" (200cm) cable
013382	Embletta Gold In-line Inductor, 4" (10cm) cable

Supplies

Natus Neuro offers a full range of neurodiagnostic supplies, featuring:

- Dedicated and knowledgeable customer service
- Streamlined order processing
- Convenient online portal (US Customers only)

Natus Neuro Store – www.natusneurostore.com

To learn more about Natus Neuro's full line of Supplies, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**



Natus Neuro

3150 Pleasant View Road
Middleton, WI 53562 USA

Tel: 1-800-356-0007

1-608-829-8500

Fax: 1-608-829-8709

www.natus.com